REBOA Catheter Convenience Set

Instructions for Use

 Prytime Medical Devices, Inc.
229 N. Main Street
Boerne, TX 78006, USA

 feedback@prytimemedical.com
www.prytimemedical.com

 US 1-210-340-0116
CAN +1-210-340-0116

REBOA CATHETER CONVENIENCE SET

CAUTION:

- USA Federal Law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).
- Prior to use, read this entire Instructions for Use.
- The below instructions, warnings, precautions, and contraindications have been incorporated from the IFUs of the constituent parts of the REBOA Catheter Convenience Set. Please consult all labeling included in the REBOA Catheter Convenience Set packaging prior to use.

DEVICE DESCRIPTION:

The REBOA Catheter Convenience Set is comprised of components compatible with the Prytime ER-REBOA™, ER-REBOA PLUS™, and pREBOA-PRO™ Catheters. The Convenience Set is utilized for gaining access, placement and fixation of the catheter.

The Convenience Set includes the following sterile items in their original packaging:

Description	Manufacturer	Device Identifier
Access Needle, 18gaX7cm, Echogenic, Sterile	Merit Medical Systems, Inc.	AD18T71W-E
Percutaneous Sheath Introducer Set, 7fr w/.035" short guidewire, Sterile	Arrow International, Inc.	CP-08703
Safety Scalpel, #11, Sterile	Aspen Surgical Products, Inc.	372611
30ml Syringe, Sterile	Becton Dickinson and Company	302832
3/4 Drape, Sterile, Disposable	Halyard Health, Inc.	47628
Cath Clamp, 5fr, Sterile	Arrow International, Inc.	CC-00005
10ml 0.9% Saline Flush, Sterile, SFR	Becton Dickenson	306553
Suture, Nylon, 2-0, Armed, Sterile	J&J Ethicon	585H

The Catheters are sold separately.

INTENDED USE:

The REBOA Catheter Convenience Set is intended to aid in gaining vascular access, placement, and securing the Catheter to the patient. The suture in this kit is intended for temporary securement of the catheter clamp to the patient's skin. The suture and clamp are to be removed immediately after the catheter is withdrawn.

CONTRAINDICATIONS:

None known.

WARNINGS:


- Read all package insert warnings, precautions, and instructions for use. Failure to do so may result in severe patient injury or death.
- The items within the Convenience Set are provided in their own Sterile, Single Use packaging. Do not reuse, reprocess or resterilize. Reuse of the devices creates a potential risk of serious injury and/or infection which may lead to death.

PRECAUTIONS:


- The Arrow Percutaneous Sheath Introducer Set contains DEHP.

INSTRUCTIONS FOR USE:**Use sterile technique.**


1. Follow Catheter Instructions for Use for placing the catheter.
2. Catheter Clamp with Fastener:
 - a. Spread wings of rubber clamp and position on catheter as required to ensure proper tip location (refer to Fig. 1).

- b. Snap rigid fastener onto catheter clamp (refer to Fig. 2).

- c. Secure catheter to patient by suturing the catheter clamp with fastener together to the skin, using the side wings to minimize the risk of catheter migration (refer to Fig. 3).

- d. Continue procedure per manufacturer's instructions.

HOW SUPPLIED

Each item within the REBOA Catheter Convenience Set is supplied sterile by the manufacturer. The individual sterile items are placed within a non-sterile over-pack. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

DEFINITIONS

Symbol	
	Read the Instructions for Use before using this product.
	Store the product appropriately in a cool, dry location.
	Do not re-sterilize this product.
	Do not reuse this product.
R_x ONLY	USA Federal Law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

Inspected Dimensions:
Folded width: 5-1/2" (14 cm)
Folded length: 7-1/2" (19 cm)

ARROW®

Catheter Clamp and Fastener Product

Safety and Efficacy Considerations:

Do not use if package has been previously opened or damaged. **Warning: Prior to use read all package insert warnings, precautions, and instructions. Failure to do so may result in severe patient injury or death.**

Do not alter the catheter clamp or fastener during use or removal.

Procedure must be performed by trained personnel well versed in anatomical landmarks, safe technique, and potential complications.

Indications for Use:

The catheter clamp and fastener are indicated for use in securing catheter to patient or as a secondary suture site.

Contraindications:

None known.

Warnings and Precautions:*

- Warning: Sterile, Single use: Do not reuse, reprocess or resterilize. Reuse of device creates a potential risk of serious injury and/or infection which may lead to death.**
- Precaution: Catheter clamp and fastener must not be attached to catheter until spring-wire guide is removed.**

A Suggested Procedure: Use sterile technique.

- Follow manufacturer's instructions for placing catheter.
- After spring-wire guide has been removed and the necessary lines have been connected or locked, spread wings of rubber clamp and position on catheter as required to ensure proper tip location (refer to Fig. 1).


Fig. 1

- Snap rigid fastener onto catheter clamp (refer to Fig. 2).


Fig. 2

- Secure catheter to patient by suturing the catheter clamp and fastener together to the skin, using side wings to minimize the risk of catheter migration (refer to Fig. 3).


Fig. 3

- Continue procedure per manufacturer's instructions.

*If you have any questions, please contact Arrow International, Inc.

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Teleflex Medical 
IDA Business and Technology Park
Athlone, Ireland

0.9% Sodium Chloride Injection, USP

BD PosiFlush™ SF Saline Flush Syringe

Description:

BD PosiFlush™ SF Saline Flush Syringes contain sterile, non-toxic, non-pyrogenic, 0.9% Sodium Chloride Injection, USP, packaged in a plastic, disposable, single use syringe that does not contain natural rubber latex or preservatives.

The package contents are sterile. Solution pH is 4.5 - 7.0. The solution osmolarity is 0.308 mOsm/mL (calc.). **Suitable for use on sterile field.**

Intended Use:

BD PosiFlush™ SF Saline Flush Syringes are intended to be used only for the flushing of indwelling vascular access devices.

Contraindications: None.

Cautions:





1. Do not use if package is not intact.
2. Do not use if syringe is damaged in any way that suggests syringe leakage.
3. Do not use if syringe tip cap is improperly affixed to or dislodged from the syringe tip.
4. Visually inspect and do not use if solution is discolored, cloudy or hazy, contains a precipitate or has any type of suspended particulate matter.
5. Do Not Re-Sterilize Before Use.
6. Verify the expiration date on the product label. Do not use if product has expired.
7. Do Not Reuse. For Single Use Only. Discard any unused portion.
8. Avoid contact between this solution and incompatible drug products. Consult appropriate compatibility literature.
9. Federal (USA) law restricts this device to sale by or on the order of a physician.
10. We guarantee the pre-filled syringe in our unopened, undamaged package to be sterile, non-toxic and non-pyrogenic.

Adverse Events: None.

Storage:

Store at 20-25°C (68-77°F); excursions permitted to 15-30°C (59-86°F); do not freeze.

Directions For Use:

1. May be used on sterile field using proper aseptic technique.
2. After removing from package, depress plunger with tip cap on to relieve the resistance between the stopper and the barrel. (Fig. 1).
3. Remove the syringe tip cap from the FLUSH syringe by twisting off. (Fig. 2).
4. Attach Blunt Plastic Cannula, if required, for flushing a needleless I.V. system with pre-slit septum. For traditional ports, attach a needle with safety engineered feature as required by the OSHA Bloodborne Pathogens Standard.
5. Hold the syringe upright and expel the air in the syringe. (Fig. 3)
Note: BD™ Blunt Plastic Cannula Utilization - When expelling air/fluid from a syringe, contents exit sideways from the cannula in two separate paths (Fig. 4). Expel air or fluid carefully, directing flow paths away from face or mucus membranes.
6. Attach the syringe to the port, valve or needleless system and flush following institution's policy and indwelling device manufacturer's recommendations.
7. Discard used syringe and any unused portion of the solution according to institution policy. DO NOT REUSE.

How Supplied:

BD PosiFlush™ SF Saline Flush Syringes are provided in sterile, peel apart packages in shelf cartons of thirty (30) unit packaged syringes, for sterile field procedures.

Reorder No. 306553 - BD PosiFlush™ SF 10 mL saline syringe with regular length plunger rod

BD PosiFlush™ Saline Flush Syringes with sterile solution and fluid path are provided in shelf containers of thirty (30) unit packaged syringes in the following fill volume / syringe size combinations.

Reorder No. 306544 - BD 3 mL syringe

Reorder No. 306545 - BD 5 mL syringe

Reorder No. 306546 - BD 10 mL syringe with regular length plunger rod

Reorder No. 306547 - BD 10 mL syringe

Reorder No. 306549 - BD 3 mL syringe packaged with separate BD Blunt Plastic Cannula

Reorder No. 306550 - BD 5 mL syringe packaged with separate BD Blunt Plastic Cannula

Reorder No. 306551 - BD 10 mL syringe with regular length plunger rod and packaged with separate BD Blunt Plastic Cannula

Becton, Dickinson and Company, 1 Becton Drive, Franklin Lakes, NJ 07417 USA

Made in Ireland bd.com

© 2016 BD. BD, the BD Logo and BD PosiFlush are trademarks of Becton, Dickinson and Company.