

# HOSPITAL REBOA KIT SUGGESTIONS

Item	Manufacturer	Model	P/N	Description
<b>Catheter</b>	Prytime	ER-REBOA™	ER7232A	7 Fr Balloon Occlusion Catheter with Blood Pressure Monitoring Arterial Line
<b>Sheath Kit</b>	Cordis	Avanti® + Sheath Introducer	402-607X	7Fr Introducer sheath requires separate 18g needle (Recommended)
<b>OR</b>				
<b>Sheath Kit</b>	Terumo	Pinnacle Precision Access System® Sheath	70-7160	7 Fr Introducer sheath and dilator, 21g micro-puncture needle, .021" nitinol guidewire
<b>Syringe</b>	Terumo	30cc Luer Lock Tip Syringe	SS-30L	30cc luer lock tip syringe
<b>Saline</b>	Hospira	0.9% Saline Solution	0409-4888-50	0.9% Saline Solution in 50mL flip-top vial
<b>Clamp</b>	Arrow®	5Fr Catheter Clamp	CC-00005	5 Fr Catheter Clamp with Fastener
<b>Scalpel</b>	Bard- Parker	No.11 Scalpel	371611	Scalpel, No. 11, Stainless Steel, Disposable, Sterile
<b>Suture</b>	Ethicon	Ethilon® 2-0 Nylon Suture	1697H	2-0 Nylon Suture 30" Reverse Cutting PSLX, 3/8C
<b>Contrast (Optional)</b>	GE	Omnipaque™ 300 mgI/ml	0407-1413-10	Contrast Solution, 10mL
<b>Needle (Recommended)</b>	Merit Medical	Merit Medical Advance® Needle	AD18T71W-EP	Merit Medical Advance® Needle, 18g, 7.0cm Echogenic needle

Item	Manufacturer	Model	P/N	Description
<b>Known Non-Compatible 7Fr Sheaths</b>	Cook	Check Flo Performer	G07431	7 Fr Introducer Sheath
	Arrow®	Percutaneous Sheath	AK-9701	7Fr Introducer Sheath
	St. Jude Medical	Fast-Cath™ Introducer With Cath-Lock & Guidewire	406702	7 Fr Introducer Sheath

Item	Manufacturer	Model	P/N	Description
<b>Known Compatible 7Fr Sheaths</b>	Merit Medical	Prelude®	PSI-7F-11-035	7FR Sheath Kit Requires separate 18g needle (Recommended)
	Merit Medical	Prelude® Pro	PRO-7F-11-035	
	St. Jude	Fast-Cath™	406108	
	Terumo	Pinnacle® R/0 II	RSB712	
	Terumo	Pinnacle® Precision w/21g micro puncture needle(SS wire)	70-7130	
	Medtronic	Input®	071101A	
	Arrow®	Super-Arrow- Flex®	CP-07711	

**Robin Cavill - West Region**  
rcavill@prytimemedical.com  
303.868.6699

**Brian Brady - South Region**  
bbrady@prytimemedical.com  
203.788.3030

**Sean Ehrhart - Great Lakes Region**  
sehrhart@prytimemedical.com  
614. 674.3345

**Seth Schreiber - Southeast Region**  
sschreiber@prytimemedical.com  
206.930.6102

**Brian Gorman - Northeastern Region**  
bgorman@prytimececal.com  
614.674.3345

**Customer Service**  
sales@prytimemedical.com  
210.340.0116